ADESIVO CIRÚRGICO GRF

1035771

BR

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EN

SURGICAL ADHESIVE GRF

BEFORE USING PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.

DESCRIPTION

SURGICAL ADHESIVE GRF (Gelatin, Resorcin, Formaldehyde) is a sterile nonabsorbable surgical adhesive composed of adhesive (Gelatin and Resorcin) and a polymerizing solution (Formaldehyde and Glutaraldehyde).

FÓRMULA

ADHESIVE:	
Gelatin	37.5%
Resorcin	12.5%
Distilled Water q.s.p.	100%
POLYMERIZING SOLUTION:	
Formaldehyde	8.3%
Glutaraldehyde	2.5%
Distilled Water q.s.p.	100%

SURGICAL ADHESIVE GRF meets all the requirements established by the Brazilian Pharmacopoeia (BR), the United States Pharmacopoeia (USP), and the European Pharmacopoeia (EP) with regards to sterilization. The gelatin used is a pharmaceutical gelatin, a proteinaceous product light yellow in color, obtained through the extraction of partially hydrolyzed collagen of animal tissue. The pharmaceutical gelatin meets the requirements of the Brazilian (BR), American (USP) and European (EP) Pharmacopoeias.

RESORCIN is a chemical compound with the molecular formula $C_7H_6O_3$, presenting a solid state and white in color. The pH (20° C) is between 4 and 6 and is soluble in water and ethanol. Resorcin may cause irritation to the skin, eyes and mucosa. Formaldehyde 40% is a chemical substance with the molecular formula CH_2O , presenting a colorless, liquid state with a pungent smell. The pH is between 2.8 and 4.0 and is soluble in water. Formaldehyde 40% is considered toxic and may cause irritation in mucosa and burns to the skin and eyes and corrosive lesions in the digestive system.

GLUTARALDEHYDE 25% is a base solution of glutaraldehyde, with the molecular formula $C_5H_8O_2$. It is a colorless liquid, with a pungent odor. Glutaraldehyde 25% may cause irritation in mucosa and skin. It is considered a toxic and corrosive substance.

INDICATIONS

SURGICAL ADHESIVE GRF is indicated for use to improve tissue hemostasis in the correction of aortic aneurism in addition to serving as suture line reinforcement. It is also indicated for use in the recovery of lacerated walls of the aorta.

ACTIONS

SURGICAL ADHESIVE GRF acts as a surgical haemostatic, used as a coadjutant in haemostatic procedures.

CONTRAINDICATIONS

SURGICAL ADHESIVE GRF is not indicated in procedures where there may be extended contact with saline solutions as the adherence of SURGICAL ADHESIVE GRF may be compromised.

WARNINGS AND CARE

The mixture of Formaldehyde and Glutaraldehyde (Polymerizing solution) is highly toxic and should not have direct contact with tissue. If direct contact occurs, immediately dry the affected area. The adhesive and Polymerizing solution mixture is not toxic.

It is recommended that gloves, bandages, instruments and other materials used during the surgical procedure remain damp to avoid SURGICAL ADHESIVE GRF adhering inadvertently to said surfaces.

Do not apply SURGICAL ADHESIVE GRF on wet tissue, as the risk of adherence may be compromised.

SURGICAL ADHESIVE GRF polymerization is rapid, therefore application is necessary soon after its components are mixed. Special care should be taken during SURGICAL ADHESIVE GRF application, since the removal of the adhesive from the tissue receiving the mixture may inadvertently cause lesions.

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